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**Stocking aid  
for compression stockings**

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1. Place your compression stocking inside the semi-circle frame and pull the top of the stocking down over the semi-circle. The heel of the stocking needs to face the back of the stocking aid.
2. Continue to slide the stocking down over the semi-circle until the heel is centered and even with the top of the semi-circle. (It may be necessary to pull the heel farther down depending on foot and stocking size.)
3. While sitting, insert your foot into stocking until your foot is completely on the floor. Helpful tip: hold the heel pocket against the frame as you slide your toes into the stocking.
4. Hold the padded handles and gently pull the stocking aid and stocking up the leg until the stocking is above the calf. Try not to overstretch the stocking.
5. Pull the stocking aid backward and down, free of the stocking.
6. Be sure the heel of the stocking is properly placed on the foot, then adjust the length and smooth out any wrinkles by stroking and smoothing with the palm of your hands.

Tigh and waist height stockings can also be donned using your stocking aid with the lower leg applied in the same manner as a knee length stocking and then pulling the upper portion into place by hand.

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